

Sampling Compatibility Codes

SLGi provides these codes to help determine whether or not analytes can be sampled on the same collection media. If the compatibility code letters match, the analytes may be sampled together. If the compatibility codes do not match, a separate tube/cassette must be used.

The applicable desorption/extraction solvents identified by compatibility codes are listed here for your information.

- A. Carbon disulfide
- B. Benzene
- C. Methanol
- D. Toluene
- E. Carbon disulfide + 5% 2-propanol
- F. Nitric acid
- G. Ethanol
- H. Carbon disulfide + 1% 2-propanol alcohol
- I. Methylene chloride + 5% methanol
- J. Carbon disulfide + 1% methanol
- K. 0.5 N Sulfuric acid
- L. 1N Sulfuric acid
- M. Carbon disulfide + 1% 2-butanol
- N. Acetonitrile/Dimethyl sulfoxide
- O. Sodium thiosulfate
- P. Carbon disulfide: XAD-2
- Q. Methanol: XAD-7
- R. Acetonitrile
- S. Acetonitrile - carbonyl free

